
DK50 2V

DK50 2V/110 DK50 2X2V/110

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IMPORTANT INFORMATION

1. CE MARKING

Products labeled with the CE mark of compliance meet the safety guidelines (93/42/EEC) of the European Union.

2. WARNINGS

2.1. General warnings

- This Installation, Operation and Maintenance Manual is a part of the appliance and must be kept with the compressor. Careful review of this manual will provide the information necessary for correct operation of the appliance.
- The safety of operating personnel and trouble-free operation of the appliance are guaranteed only if original parts are used. Only accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer can be used.
- If any other accessories or consumable materials are used, the manufacturer cannot be held responsible for the safe operation of the appliance. This guarantee does not cover damages originating from the use of accessories or consumable material other than those specified or suggested by the manufacturer.
- The manufacturer guarantees the safety, reliability and function of the appliance only if:
 - Installation, new settings, amendments, extensions and repairs are performed by the manufacturer or its representative, or a service provider authorized by the manufacturer
 - The appliance is used in accordance with this Installation, Operation and Maintenance Manual
- The manufacturer reserves all rights for the protection of its wiring diagrams, methods and names.
- Translation of Manual for Installation, Operation and Maintenance is carried out in accordance with the best knowledge. In the case of ambiguities, the Slovak version of the text prevails.

2.2. General safety warnings

The manufacturer developed and designed the equipment in such a way so that any risks were excluded if it is used according to intention. The manufacturer considers it to be its obligation to describe the following safety measures in order to exclude residual damages.

- Operation of the appliance must be in compliance with all local codes and regulations.
- Original packaging should be kept for the return of the appliance. Only the original packaging ensures protection of the appliance during transport. If it is necessary to return the appliance during the guarantee period, the manufacturer is not liable for damages caused by improper packaging.
- Each time the appliance is used, the operator must make sure that it is functioning correctly and safely.
- The user must fully understand the operation of the appliance.
- The product is not intended for operation in areas with a risk of explosion.
- If any problem occurs during use of the appliance, the user must inform his supplier immediately.

2.3. Electrical system safety warnings

- The appliance must be connected to earth (grounded).
- Before the appliance is plugged in, make sure that the mains voltage and mains frequency stated on the appliance are the same as the power mains.
- Prior to putting into operation it is necessary to check for possible damage of the equipment and connected air and electric distributions. Damaged pneumatic and electric lines must be immediately replaced.
- Immediately disconnect the appliance from the mains (pull out mains plug) if a technical failure occurs.
- During repairs and maintenance, ensure that:
 - The mains plug is pulled out from the socket
 - Pressure pipes are vented and pressure is released from the air tank.
- The appliance must be installed by an approved, qualified technician.

3. ALERT NOTICES AND SYMBOLS

In the Installation, Operation and Maintenance Manual and on packaging and product, the following labels or symbols are used for important information:

	Information, instructions and cautions for the prevention of damage to health or materials
	Caution! Dangerous electrical voltage
	Read the user manual!
	CE mark of compliance
	Compressor is remote-controlled and may start without warning
	Caution! Hot surface
	Earth (ground) connection
	Terminal for ground connection
	Fuse
	Alternating current
	Handling mark on package – Fragile, handle with care
	Handling mark on package – This way up (vertical position of cargo)
	Handling mark on package – Protect against moisture
	Handling mark on package – Temperature during storage and transport
	Handling mark on package – Limited stacking
	Mark on package – Recyclable material

4. STORAGE AND TRANSPORT

The compressor is shipped in cardboard that protects the appliance from damage during transport.

Caution! For transport, always use the original packaging and secure the compressor in the upright position.

Protect the compressor from humidity and extreme temperatures during transport and storage. A compressor in its original packaging can be stored in a warm, dry and dust-free area. Do not store near any chemical substances.

Keep packaging material if possible. If not, please dispose of the packaging material in an environmentally friendly way and recycle if possible.

Caution! Before moving or transporting the compressor, release all the air pressure from the tank and hoses and drain the condensed water.

5. TECHNICAL DATA

	DK50 2V	DK50 2VS	DK50 2V/110	DK50 2V/110S	DK50 2x2V/110	DK50 2x2V/110S
Nominal voltage / (*) frequency V / Hz	230 / 50 230 / 60 3x400/50	230 / 50 230 / 60 3x400/50	230 / 50 230 / 60 3x400/50	230 / 50 230 / 60 3x400/50	230 / 50 230 / 60 3x400/50	230 / 50 230 / 60 3x400/50
Efficiency of compressor at over- pressure 5 bar Lit.min ⁻¹	140	140	140	140	280	280
Efficiency of compressor with dryer at over-pressure 5 bar Lit.min ⁻¹	140	140	105 (**)	105 (**)	220(**)	220 (**)
Efficiency of compressor with KJF at over-pressure 5 bar Lit.min ⁻¹	140	140	140	140	280	280
Maximal current A	7,5 8,5 3,6	7,7 8,7 3,8	7,5 8,5 3,6	8 9 4,1	15 17 7,2	15,5 17,5 7,7
Maximal current of compressor with dryer A	7,8 8,8 3,9	8 9 4,1	7,8 8,8 3,9	8,3 9,3 4,4	15,3 17,3 7,5	15,8 17,8 8
Motor performance kW	1,1	1,1	1,1	1,1	2x1,1	2x1,1
Air tank capacity Lit.	25	25	110	110	110	110
Pressure range bar	5,0 – 7,0	5,0 – 7,0	5,0 – 7,0	5,0 – 7,0	5,0 – 7,0	5,0 – 7,0
Maximum operating pressure of safety valve bar	8,0	8,0	8,0	8,0	8,0	8,0
Sound level L _{pFA} [dB]	71	53	71	52	73	53
Mode of operation of compressor or compressor with s KJF	continual S 1	continual S 1	continual S 1	continual S 1	continual S 1	continual S 1
Mode of operation of compressor with dryer MONZUN - M1a with dryer MONZUN – M2	intermittent S 3-60%	intermittent S 3-60%	continual S 1	continual S 1	continual S 1	continual S 1
Condensation unit	KJF-1	KJF-1	KJF-1 / KJF-2	KJF-1 / KJF-2	KJF-2	KJF-2
Dimensions of compressor / of compressor with dryer w x l x h mm	460x460x690/ 460x535x690	560x665x860	1090x425x815/ 1090x605x815	1200 x 720 x 990	1090x425x815/ 1090x605x815	1200 x 720 x 990
Dimensions of compressor packaged in cardboard in box w x l x h mm	550x580x760	600x680x895	1180x720 x1030	1330x840 x1185	1180x720 x1030	1330x840 x1185
Weight of compressor / of compressor with dryer kg	52 /62	88 / 99	70 / 90	145 / 165	98 / 118	173 / 193
Weight of compressor / of compressor with dryer Packaget in cardboard in box kg	57 / 67	94 / 105	91 / 110	227 / 246	120 / 142	239 / 257
Drying point of compressor Atmospheric condensation point	-20°C	-20°C	-20°C	-20°C	-20°C	-20°C
Version EN 60 601-1	Appliance of type B, class I					

(*) When ordering, please specify the version of the compressor

(**) Part of the efficiency of compressor (cca-20%) is used for regeneration of the dryer

- Weight of compressor with KJF1 increases about 3kg, with KJF2 increases about 4kg

- Air outgoing from KJF or dryer M1and / M2 is filtered using 5µm filter

Climatic conditions during storage and transport
Temperature : -25°C to +55°C, 24 h to +70°C
Relative air humidity : 10% to 90 % (no condensation)

Climatic operation conditions
Temperature : +5°C to +40°C
Relative air humidity : 70%

6. PRODUCT DESCRIPTION

6.1. Model variations and their uses

Compressors are the source of clean, oil-free compressed air used to drive dental appliances and equipment.

Compressors models are designed for the following uses:

Dental compressors DK50 2V - are designed for independent placement of the compressor in any area.

Dental compressors DK50 2V/K - are designed for independent placement of the compressor in any area and feature a condensation and filtration unit (KJF1).

Dental compressors DK50 2V/M - are designed for independent placement of the compressor in any area and feature an adsorption air dryer M1a.

Dental compressors DK50 2VS - feature soundproof boxes suitable for placing in the dentist's surgery.

Dental compressors DK50 2VS/K - feature soundproof boxes with a condensation and filtration unit (KJF1).

Dental compressors DK50 2VS/M - feature soundproof boxes and an adsorption air dryer M1a.

Dental compressors DK50 2V/110 and DK50 2x2V/110 - These compressors can be situated in areas where their operation and noise will not disturb their immediate surroundings. They are capable of providing compressed air for several dentists' surgeries.

Dental compressors DK50 2V/110/K and DK50 2x2V/110/K - are equipped with a condensation and filtration unit.

Dental compressors DK50 2V/110/M and DK50 2x2V/110/M - are equipped with an adsorption air dryer M2.

Dental compressors DK50 2V/110S and DK50 2x2V/110S - are placed in compact soundproof boxes and are capable of providing compressed air for several dental units.

Dental compressors DK50 2V/110S/K and DK50 2x2V/110S/K - are placed in compact soundproof boxes and are equipped with a condensation and filtration unit.

Dental compressors DK50 2V/110S/M and DK50 2x2V/110S/M - are placed in compact soundproof boxes and are equipped with an adsorption air dryer M2.

Box S110 - serves for the decrease of level of compressor noise. Case with compressor may be placed directly in the dental consulting room or laboratory as the part of dental furniture.


DK50 2V

DK50 2VS

DK50 2x2V/110

DK50 2V/110

S110
 DK50 2V/110S
 DK50 2x2V/110S

KJF1

KJF2

MONZUN M1a

MONZUN M2

Without additional filtration equipment, the compressed air from a compressor is not suitable for the operation of breathing appliances or similar equipment.

7. FUNCTION

Compressor (Fig.1, Fig.2)

The compressor (1) draws in air through a filter (8) and compresses it through a check valve (3) into an air tank (2). The connected apparatus draws the compressed air from the air tank until the pressure drops to a default preset level on the air-pressure switch (4) switching the compressor on. The compressor again compresses air into the nozzle until the maximum pressure is reached and the compressor switches off. After compressor aggregate is switched off, pressure hose shall be pressure-release solenoid valve (13). Safety valve (5) prevents the pressure in air chamber from rising above the maximal allowed value. The drain valve (7) releases the condensate from the air nozzle. Compressed, clean air free from oil traces is stored in the air tank ready for use.

Compressor with dryer M1a (Fig.3)

The compressor (1) draws in air through a filter (8) and sends it to the air dryer in compressed form. The air proceeds through the cooler (15) and through the dryer chamber (9) with desiccant (16). The moisture is captured in an output filter (17) and check valve (3), sending clean, dry air into the air tank (2). Adsorber shall be regenerated after every switching off of compressor using pressure switch during dryer chamber depressurisation. Air escapes from adsorption chamber via open solenoid valve (14), while it is concurrently blown with dried air. The level of drying of the passing-through air is monitored by a humidity sensor (23), that shall open the regeneration valve (24) in the case of registering humidity exceeding the set humidity and at the time of a pause of the compressor it shall automatically regenerate the drying agent in the chamber. The entrapped water from the chamber shall be pushed out via a valve. Dry, clean compressed air free from oil traces is stored in the air tank ready for use.

Compressor with dryer M2 (Fig.4)

The compressor unit (1) pulls air in via the intake filter (8) and supplies compressed air to the cooler (15). The compressed air then passes to the condensate separator (43) through the drying chamber (9) with the absorber (16) where moisture is removed and then out via the outlet filter (17) and the check valve (3) and into the air tank (2). The absorber regenerates when the drying chamber is being evacuated, always after the compressor is switched off by the pressure switch. The air is then released from the absorption chamber via the open solenoid valve (14) when it is cleared using dry compressed air. Air drying occurs in one chamber while regeneration occurs in the other chamber. The mode of each chamber changes in regular cycles and the drying and regeneration processes are then performed in the alternate chamber. Compressed, dry and clean air without any trace of oil is then ready for additional use in the air tank.

Compressor with condensation and filtration unit KJF1, KJF2 (Fig.5)

The compressor (1) draws in air through a filter (8) and compresses it through a check valve (3) into an air tank (2). The compressed air from the nozzle flows through a cooler (10) that cools the compressed air. The condensed moisture is trapped in the filter (11) and automatically separates as condensate (12). Dried, clean compressed air, free from oil traces, is ready for use.

Compressor box (Fig.6, Fig.7)

The box is soundproof and allows sufficient exchange of cooling air. It can be placed in a dentist's office. Fan under compressor aggregate provides for the compressor cooling. It is in operation at one time with compressor engine, or when temperature in the casing rises above 40°C. After the space in the casing cools down under ca 32°C, fans shall get automatically switched off. Door of the casing with right opening may be changed to left opening (see Chapter 9).


Make sure that nothing impedes the free flow of air under and around the compressor. Never cover the hot air outlet on the top back side of the case.


If placing the compressor on a soft floor such as carpet, create space for ventilation between the base and floor or the box and floor, e.g. underpin the footings with hard pads.

Fig.1 – Compressor

1. Compressor motor
2. Air tank
3. Check valve
4. Pressure switch
5. Safety valve
6. Manometer
7. Drain valve
8. Input filter
9. Dryer chamber
10. Pipe cooler
11. Output filter
12. Condenser outlet
13. Solenoid valve
14. Output solenoid valve of dryer
15. Cooler
16. Dessicant
17. Output filter
18. Sieve
19. Box fan
20. Stopper
21. Control valve
22. Bottle
23. Humidity sensor
24. Regeneration valve
25. Box
26. Lock
27. Connecting reinforcement
28. Wall stopper
29. Switch
30. Manometer
31. Magnetic bottle holder
32. Door hinge
33. Wheels
34. Socket on the box
35. Silencer
36. Hose PUR ø8 /ø6
37. Screwing with a tap
38. Straight screwing
39. Power supply cord
40. Hose of manometer
41. Compressor Fan

Fig.2 – Compressor

22. Bottle
23. Humidity sensor
24. Regeneration valve
25. Box
26. Lock
27. Connecting reinforcement
28. Wall stopper
29. Switch
30. Manometer
31. Magnetic bottle holder
32. Door hinge
33. Wheels
34. Socket on the box
35. Silencer
36. Hose PUR ø8 /ø6
37. Screwing with a tap
38. Straight screwing
39. Power supply cord
40. Hose of manometer
41. Compressor Fan
42. Angular screwing ¼"M-8/6"
43. Water separator
44. Rectification screw
45. Door pin
46. Compressor handle
47. Lock S110
48. Screw M5
49. Box gas springs
50. Drain valve
51. Fitting of output pressure hose
52. Angular screwing 3/8" MF
53. Handle S110
54. Vessel
55. Safety valve
56. Magnetic bottle holder S110
57. Vessel S110

Fig. 3 - Compressor with air dryer MONZUN – M1a


Fig. 4 - Compressor with air dryer MONZUN – M2


Fig.5 - Compressor with condensation and filtration unit KJF


Fig.6 - Box





Fig.7 - Box S110

INSTALLATION

8. USE

- The appliance must be installed and operated in a dry, well ventilated and dust-free area where ambient temperature is within the range of +5°C to +40°C and relative air humidity does not exceed 70%. Otherwise, failure-free operation of the compressor cannot be guaranteed. The compressor must be installed so that it is accessible at all times for operation and maintenance. Please ensure that the appliance label is accessible.
- The appliance must stand on a flat, sufficiently stable base. See paragraph 5 (Technical data) when positioning or lifting the compressor.
- Compressors cannot be exposed to outdoor environments. The appliance cannot be used in moist or wet environments. Do not use the compressor in the presence of explosive gases, dust or combustible liquids.
- Before connecting the compressor to medical equipment, the supplier must confirm that it meets all requirements for its use. Refer to the technical data of the product for this purpose. When a unit is to be built-in, classification and evaluation of compatibility must be done by the manufacturer or supplier of the product to be used.
- Any use other than that described in this manual is not covered by the guarantee, and the manufacturer is not liable for any damages that may result. The operator/user assumes all risk.

9. INSTALLATION

Only qualified personnel can install and start up the appliance and train operating personnel in its correct use and maintenance. Installation and training of all operators shall be confirmed by the installer's signature on the certificate of installation.

Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.

Caution! When in operation, the compressor is hot. Burns or fire may result if contact is made by the operator or any flammable material.

Electric cord for connection to electric mains and air hoses may not be broken. The power cord may not be exposed to pulling, pressure and excessive heat.

9.1. Placement of the compressor

Handling


Fig.8 Unpacking


Fig.9 Unpacking

Dental compressor DK50 2V, DK50 2V/M, DK50 2V/110, DK50 2V/110/M, DK50 2x2V/110, DK50 2x2V/110/M (Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Connect output hose with end-piece to the appliance. Plug the mains cord plug into a socket. Place the drain hose into a container prepared in advance.

Dental compressor in box DK50 2VS (Fig.6, Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Place the wall-mounted stopper (28) - 2 pcs onto the compressor housing in the rear top part of the housing and put the housing onto a required place. The stoppers provide a sufficient distance of the housing from a wall for thorough ventilation. For setting up the compressor in housing you must open the door on the housing using the attached key and remove connecting reinforcement (27) in the front bottom part of housing. If necessary, the door may be disassembled using door hinge (32). Connect the compressor via distributions prepared in advance in a floor according to the installation plan or via holes in the rear part of the housing (Fig.10). Protrude pressure hose via a hole in housing and connect it to the appliance in a suitable way. Grasp the compressor at its handle and put it into the casing using built-in wheels (33). Fix the hose (40) of a manometer (30) in a hose into the fast-on coupling on a compressor, put the connecting reinforcement (27) back and connect the pressure hose to a compressor. Insert the electric power cord of a compressor into a socket (34) on a housing. By slight rotation of rectification screws (44) set the correct position of door against the casing frame. When closing the door the pin (45) on the door must easily snap in the opening in the casing frame. Close the housing doors and duly lock the lock (26). Connect the mains plug into the mains socket.

It is not allowed to leave the key in a lock! It must be saved against non-instructed persons!

Dental compressor in box DK50 2VS/M (Fig.6, Fig.8)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Place the compressor into housing similarly as in the previous paragraph. Prior to placing the compressor into housing, protrude house for condensate drain via hole in housing and connect it to a bottle (22). Magnetic holder (31) with a vessel (22), for entrapping the condensate from the dryer may be fitted onto any vertical part of casing, or from front on its door. When fixing the holder with a vessel at the housing side it is necessary to consider a space of at least 11 cm between the housing and furniture. Distance smaller than the specified one may cause problem with handling of the vessel.

Dental compressor in box DK50 2V/110S, DK50 2x2V/110S(Fig.7, Fig.9)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Place the wall-mounted stopper (28) - 2 pcs onto the compressor housing in the rear top part of the housing and put the housing onto a required place. The stoppers provide a sufficient distance of the housing from a wall for thorough ventilation. Dismount the door, fixed using 4 screws M5 (48) located in the tapered edges of the door and disconnect earthing wire. If necessary, the top panel of the housing may be opened by releasing the locks (47) by moderate rotation of screw driver according to pictogram and lifting using a handrail (53). Gas springs (49) hold it in an open position. Insert the compressor into the housing so that the compressed air outlet was oriented towards the operators and **so that a gap of at minimally 50 mm was between electric motor and plastic foam on the rear panel of the housing.** Orientate the output pressure hose to the rear part of the compressor. Connect the power cord of the compressor (39) into the socket (34) on the electric panel of the housing. Place the threaded joint with a tap (37) into the hole on

the side of the housing and place PUR hose with $\text{Ø}8 / \text{Ø}6$ (36). Blind the hole on the opposite side of the housing with a plug with $\text{Ø}15.5$. (*Selection of the side for placing threaded joint with a tap depends upon the decision of a customer*). The other end of the hose (36) place into the straight threaded joint (38) on the air chamber. Place a hose (40) leading from manometer into the angular threaded joint on the air chamber (42). Protrude the pressure hose via a hole in the housing and connect it to the appliance in a suitable way. Connect the power cord plug leading from the housing into the mains socket. Insert the connector of earthing wire onto the door and screw it onto the housing. Cover the holes after the screws using white end caps with $\text{Ø}11$. Close the top panel of the housing and lock the locks.

Dental compressor in box DK50 2V/110S/M, DK50 2x2V/110S/M (Fig.7, Fig.9)

After removing all packaging material, place the product on the floor and remove stabilization parts X and Y (Detail A). Place the wall-mounted stopper (28) - 2 pcs onto the compressor housing in the rear top part of the housing and put the housing onto a required place. The stoppers provide a sufficient distance of the housing from a wall for thorough ventilation. Dismount the door, fixed using 4 screws M5 (48) located in the tapered edges of the door and disconnect earthing wire. If necessary, the top panel of the housing may be opened by releasing the locks (47) by moderate rotation of screw driver according to pictogram and lifting using a handrail (53). Gas springs (49) hold it in an open position. Prior to placing the compressor into housing, protrude house for condensate drain via hole in housing and connect it to a bottle (57). Magnetic holder (56) with a vessel (57), for entrapping the condensate from the dryer may be fitted onto any vertical part of casing, or from front on its door. When fixing the holder with a vessel at the housing side it is necessary to consider a space of at least 16 cm between the housing and furniture. Distance smaller than the specified one may cause problem with handling of the vessel. Insert the compressor into the housing so that the compressed air outlet was oriented towards operators and **so that ventilator of a dryer was inserted into the venting tunnel in a housing**. Orientate the output pressure hose to the rear part of the compressor. Connect the power supply cord of compressor (39) to the socket (34) on electric panel of the housing. Blind the holes at the sides of the housing using a plug with $\text{Ø}15.5$. Place a hose (40) leading from manometer into the angular threaded joint on the air chamber (42). Protrude the pressure hose via a hole in the housing and connect it to the appliance in a suitable way. Connect the power cord plug leading from the housing into the mains socket. Put on the connector of earthing wire onto the door and screw it onto the housing. Cover the holes after the screws using white end caps with $\text{Ø}11$. Close the top panel of the housing and lock the locks.

Dental compressor DK50 2V/110/M, DK50 2x2V/110/M, DK50 2V/110, DK50 2x2V/110 which was not delivered together with the case S110 (Fig.7, Fig.9)

Prior to the placing of compressor into the case, it is necessary to make the following works at the compressor:


- Remove tap for the drainage of condensate (15) from air chamber (bonded joint) and to replace it with straight fitting $\text{Ø}8 - \frac{1}{4}$ " (12) (bond). – for version of compressor without dryer
- Remove the fitting (16) for the connection of output pressure hose (bonded joint) and place it into angular fitting (17) $\frac{3}{8}$ " M/F (bond). Then fix this unit to the original place of fitting (bond) so that the output of air was directed downwards.
- Remove the plug (bonded joint) and to replace it with angular fitting (13) $\frac{1}{4}$ M-8/6 (bond).

Place the compressor into the housing similarly as in the previous paragraphs for compressor in a housing DK50 2V/110S/M, DK50 2x2V/110S/M or compressor in a housing DK50 2V/110S, DK50 2x2V/110S.

While closing the upper board, it is necessary to increase the care – the risk of squeezing the fingers of a hand.

After closing the lid of the case it is necessary to secure the fast-on elements!


- 1 - Outline of the case
 - 2 - Outline of the base
 - 3 - Front part of the case - door
 - 4 - Pressurised air connector G3/8"
 - 5 - Supply voltage
 - 230V/50(60)Hz 3Gx1.5x1000
 - 400V/50(60)Hz 5Gx1.0x3000
- all dimensions in millimeters
- wall distance min. 100mm

9.2. Compressed air outlet

(Fig.11)

Lead the pressure hose from the output of compressed air (1) to the appliance – dental set.

9.3. Electrical connection

Fig.11

Plug the electrical cord into the mains.

The appliance is equipped with a grounded plug. Make sure this connection complies with local electrical codes. The mains voltage and frequency must comply with the data stated on the appliance label.

(Fig.12)

- Keep the socket easily accessible to ensure that in an emergency the appliance can be safely disconnected from the mains.
- Connection to the power distribution box must be max.16 A.
- The connection of the earth ground pin \varnothing 6mm (1) with other appliances must be completed in accordance with local electrical codes. The female socket (2), which is not included in the standard set, is an optional accessory.

Fig.12


Electrical cable may not contact the hot parts of a compressor. Insulation could be damaged!



If any electrical cord or air hose is damaged it must be replaced immediately.

9.4. Change in the door opening




- Disassemble the door, rectification screw (3) and the holder (2) of the hinge D (4).
- Mount holder of the hinge D to the left side of the casing.
- Rotate the door by 180°.
- Insert a spacer (5) between the hinge H(1) and the bottom side of the door
- Mount the door.
- Disassemble the lock (6) on the door, rotate it by 180°.
- Disassemble the latch (7) on the lock, rotate it by 180°.
- Mount the lock to the door.

B


C


10. WIRING DIAGRAMS

1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE


- M1 Motor of compressor
- EV1 Fan of compressor
- EV3* Fan of dryer
- EV4 Fan of box
- FA1 Breaker
- ST1, ST2 Thermo switch
- SP Pressure switch
- YV1 Solenoid valve of compressor
- YV2* Solenoid valve of dryer
- X1, X2 Terminal
- X3 Terminal with fuses
- XC Connector
- SRH* Humidity sensor
- HL1 Glowlamp
- Q Switch
- CB1 Capacitor

DK50 2V
DK50 2V/110


1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE

- SP Pressure switch
- EV1, EV2 Fan of compressor
- ST1 Thermo switch
- EV3* Fan of dryer
- M1, M2 Motor of compressor
- YV1 Solenoid valve of compressor
- FA1, FA2 Breaker
- YV2*, YV4* Solenoid valve of dryer - OUT
- KM1, KM2 Contactor
- YV3* Control valve
- X1 Terminal with fuses
- KT Time relay
- PS 2V2 PCB


DK50 2V/110

1/N/PE ~ 230 V 50/60 Hz
ELECTRIC OBJECT OF 1 ST CAT
TYPE B


DK50 2V / 110 / M

1/N/PE ~ 230 V 50..60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE


DK50 2x2V/110

1/N/PE - 230 V 50/60 Hz
ELECTRIC OBJECT OF 1 ST CAT
TYPE B

DK50 2 x 2V / 110 / M


3/N/PE - 400/230 V 50 Hz
MAINS TN-S [TN-C-S]
ELECTRIC OBJECT OF 1st CAT.
B TYPE

DK50 2x2V/110

- M1, M2 Motor of compressor
- KM1, KM2 Contactor
- FA1, FA2 Breaker
- ST1 Thermo switch
- SP Pressure switch
- PS 2V2 PCB
- EV1, EV2 Fan of compressor
- EV3* Fan of dryer
- YV2*, YV4* Solenoid valve of dryer - OUT
- YV3* Control valve
- KT Time relay
- X1 Terminal with fuses

3/N/PE - 400/230 V 50.60 Hz
MAINS TN-S [TN-C-S]
ELECTRIC OBJECT OF 1 ST CAT
TYPE B


DK50 2x2V/110/M

DK50 2V DK50 2V/110


3/N/PE ~ 400/230 V 50 Hz
MAINS TN-S [TN-C-S]
ELECTRIC OBJECT OF 1st CAT.
B TYPE

DK50 2V , * DK50 2V/M


- Q Switch
- HA, HL1 Glow lamp
- KM1, KM2, KM3 Contactor
- EV4, EV5 Case fan
- X1, X2 Fuses terminal
- XC1 Connector
- ST1, ST2 Thermal switch
- SP Pressure switch
- EV1 Fan of compressor
- EV3* Fan of dryer
- YV1 Solenoid valve of compressor
- YV2* Solenoid valve of dryer
- SRH* Humidity sensor
- FA1 Breaker
- M1 Motor

3/N/PE ~ 400/230 V 50, 60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE


DK50 2VS

3/N/PE ~ 400/230 V 50, 60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE

S110

1/N/PE ~ 230 V 50, 60 Hz
ELECTRIC OBJECT OF 1st CAT.
B TYPE

S110

11. FIRST OPERATION

(Fig.13)

- Make sure that all stabilizers used during transport were removed.
- Check that all pressurized air line connections are secure.
- Connect to the mains.
- Start compressor at pressure switch (2) by turning switch (3) to position "I."
- For kompressor in the box turn the switch (29) (Fig.6, Fig.7) at the front part of the soundproof box to the position "I" – green light indicates that the appliance is on.

Compressor - At first operation the air tank is pressurized until it reaches a preset level when the compressor automatically switches off. As the air is used, the compressor works in automatic mode, switched on or off by the pressure switch.

Compressor with dryer - In addition, adsorption dryer takes out humidity from the previous compressed air during operation in the device and it blows out the entrapped condensate via condensate outlet on the dryer, which is audible as a short hiss when compressor stops, or during the operation while switching the dryer chambers (for M2).

Compressor with condensation and filtration unit - Model KJF filters and dehumidifies the air and automatically releases condensed liquid through the filter's discharge valve.

The compressor is not equipped with an emergency power supply.

OPERATION

In case of emergency, disconnect the compressor from the mains (pull out the mains plug).

**The compressor has hot surfaces.
Burns or fire may result if contact is made.**

During prolonged operation of the compressor, the temperature in the box may increase to over 40°C. At this point the cooling fan automatically switches on. After cooling the space to under 32°C, the ventilator switches off.

Automatic start: when pressure in the tank drops to the pressure switch's lower limit level, the compressor automatically switches on. The compressor automatically switches off after reaching the pressure switch's upper limit level.

Compressor with dryer M1a

A correct function of the drier depends on the compressor's operation and no attendance is required. The pressure vessel need not be sludged, because the pressure air entering the air chamber is already dried.

- For correct operation of a drier it is necessary to:
- Observe the operation of the compressor in intermittent mode up to 60%. However, the time of continuous operation of the compressor should not be longer than 10 minutes.
- It is forbidden to alter the working pressures of pressure switch set by manufacturer. The operation of the compressor at working pressure lower than the switching pressure demonstrates the overload of the compressor (high air consumption) by the appliance, leakages in pneumatic distributions, failure of aggregate or drier.
- Leave the compressor connected to electric mains (do not switch off pressure switch and do not disconnect power cord) – if the humidity sensor assesses air as "insufficiently dry", it shall regenerate the drying chamber refill by intake of air from air chamber in which pressure drops and it may come to switching on the compressor several times. After the drier is regenerated, the operation of the compressor shall be automatically halted.
- Should regeneration takes longer than 1 hour and there is no air consumption, it shall be necessary to verify whether:
 - Regeneration takes place,
 - Air goes out from the outlet of solenoid valve of the drier via condensate drain,
 - There is no failure of the compressor or drier

Compressor with dryer M2

- It is forbidden to alter the working pressures of pressure switch set by manufacturer. The operation of the compressor at working pressure lower than the switching pressure demonstrates the overload of the compressor (high air consumption) by the appliance, leakages in pneumatic distributions, failure of aggregate or drier.
- Prior connecting drier to air chamber, that was used with compressor without drier, it is necessary to clean interior surface of air chamber and perfectly remove condensed liquid. Then interconnect electric part of drier with compressor according to wiring diagram in accord with valid regional regulations.

12. SWITCHING THE COMPRESSOR ON

(Fig.13)

Switch on the compressor at the pressure switch (2) by turning the knob (3) to position “1.” (for compressor in the box switch (29) Fig.6, Fig.7), on the front part of the compressor box), The compressor sends pressurized air to the air tank. As the compressed air is used, the pressure in the air nozzle drops to a preset level, the compressor switches on and the air nozzle fills with compressed air. After reaching the cutoff pressure the compressor turns off automatically and the cycle is repeated. Check the value of switching-on and switching-off pressure on pressure gauge. The values may be within a tolerance of ±10%. Air pressure in air chamber must not exceed maximal permitted operation pressure.


Fig.13

Never tamper with the pressure switch (2). Adjustments are not allowed. The pressure switch (2) has been set by the manufacturer and further setting of switching on and off pressure may be carried out only by a qualified expert trained by the manufacturer.

After putting the compressor with a drier M1a to operation (after installation or several-day break) the regeneration of the drier may be started. The compressor is put into operation automatically also without pressure air take-off by the appliance, and thus increase pressure in the air jet into cut-off pressure when it turns off. Then the dryer regeneration takes place (passing of the air from air jet through dryer chamber). The air jet pressure decreases to closing pressure, the compressor turns on and increases pressure into air jet to cut-off pressure and turns off. The process of turning on and off of the compressor repeats in this way until the dryer is not regenerated sufficiently. This level is controlled in the dryer through the built-in humidity monitor – hygostat. The regeneration process can take several minutes (5 – 15 min) – for the new or previously already regenerated dryer during the previous operation, or several tens of minutes (30 – 120 min) – for the dryer being previously in its operation „overloaded“ with water vapour (for example due to the activity regime of the compressor outside the permitted extent, during the work in the environment with high relative humidity, etc.). After finishing regeneration the whole process is stopped automatically.

Note!


The mentioned values of the regeneration times are valid for compressor activity regime at closed outlet valve, which means without pressure air take-off by the appliance. In case of air take-off the mentioned times are prolonged.

In the case the compressor does not terminate the automatic regeneration of a drier within 120 min, it shall be necessary to turn to your supplier or servicing centre.

The table 1 states the time of regeneration cycle of the compressor with a drier M1a (under the condition that no air is taken from the compressor by an appliance).

(Table 1)

Compressor	Jet size	Closing pressure - cut-off pressure	Compressor activity time	Compressor break time - dryer regeneration
DK50 2V/M1a, DK50 2V S/M1a	Ø 0,7 mm	5 – 7 bar	about 30 – 40 s	about 90 – 110 s

Used jet :
DK50 2V- 0,7 mm

In case of deviation from described activity please verify:

Connection of the solenoid valve – compare arrow orientation on the valve body, which also marks correct air flow direction when regenerating

Possible leakage via compressor untightness – close the outlet valve on the compressor, turn on the compressor and leave it in operation until it turns off at cut-off pressure. Turn off the circuit breaker and watch pressure on manometer, which means the pressure in air jet. The pressure decrease can not be bigger than 0.2 bar in 2 hours.

The use of the correct jet – a jet is mounted between solenoid valve and opening (see Table 1). When assembling it is important to keep jet orientation – aim the surface with a screwed cone at the valve.

Size of regeneration nozzle for dryer M2
 (Table 2)

Compressor	Jet size	Closing pressure - cut-off pressure	Amount of dried air	Amount of regenerated air
DK50 2V/110/M2	Ø 0,7 mm	5 – 7 bar	140 l/min.	25 l/min.
DK50 2x2V/110/M2	Ø 1 mm	5 – 7 bar	280 l/min.	50 l/min.

MAINTENANCE

13. MAINTENANCE SCHEDULE

Notice!

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

Maintenance that must be performed	Chapter	Time interval	Performed by
<ul style="list-style-type: none"> • Release condensate <li style="padding-left: 20px;">Compressor without air drier <li style="padding-left: 20px;">At high air humidity <li style="padding-left: 20px;">Compressors with air drier <li style="padding-left: 20px;">Compressors with condensation unit : <ul style="list-style-type: none"> - from filter - from pressure vessel 	14.1	1 x week 1 x day 1 x week, check function 1 x week, check function 1 x week	operating staff operating staff operating staff operating staff operating staff
• Check safety valve	14.2	1 x year	qualified technician
• Replacement of the input filter and prefilter	14.3	1 x 2 years or after 5000 hours	qualified technician
• Replacement of filter in dryer	14.4	1 x year	qualified technician
• Replacement of filter in condensation unit	14.5	1 x year	qualified technician
• Replacement of the buoy in the water separator	14.6	1 x year	qualified technician
<ul style="list-style-type: none"> • Check tightness of joints • Overall examination of device 	Service documentation	1 x year	qualified technician
• Perform "Repeated Test" according to EN 62353	13	1 x 2 years	qualified technician

14. MAINTENANCE

Repair work beyond normal maintenance can be performed only by qualified personnel or the manufacturer's representative.

Use only spareparts and accessories approved by the manufacturer.

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).


For permanently high efficiency of drying, it is necessary to maintain the whole appliance, and mainly ventilator clean – regularly clean the surface of ventilator and cooling fins of cooler.

TO ENSURE THAT THE COMPRESSOR WORKS CORRECTLY, PERFORM THE FOLLOWING MAINTENANCE TASKS AT REGULAR INTERVALS (CHAPTER 13).:

14.1. Condensation drain valve

Compressors (Fig.14, Fig.15)

During regular use, release condensation from the pressure tank. Switch off the compressor at the mains. Reduce air pressure in the appliance to max. 1 bar by releasing air via a connected device. Place the hose with the drain valve into a container prepared in advance (for compressors DK50 2V/110 and DK50 2x2V/110 orientate the vessel towards the drain valve, for compressors DK50 2V place the vessel under release valve) and open the drain valve (1). Wait until condensation is fully drained from the pressure tank. Close drain valve (1).

Compressors with condensation and filtration unit (Fig.18)

During regular use, condensation is automatically released via the release valve of the condensation unit filter. To check that the automatic drain is working properly, open the valve (4) of the drain vessel (2) by turning to the left. Release a small amount of condensate from the vessel. Close the valve (4) by turning to the right.

Fig.15
DK50 2V

Compressors with air dryer

In the case of a regular operation condensate is automatically excreted via air dryer and it is entrapped in a bottle. Take out the bottle from a holder, release blowing-out suppressor and pour out the condensate. If necessary, it is possible to connect the set for condensate discharge onto the condensate outlet (see Chap. PARTS LIST - Auxiliary Equipment).

For versions of a compressor with a housing it is necessary to open the housing prior to the following checks.

For DK50 2VS - unlock the lock on the door and open the door of the housing (Fig.6)

For DK50 2V/110S, DK50 2x2V/110S unlock the locks on the top panel of the housing and lift (Fig.7)

14.2. Safety valve check

(Fig.13)

When the compressor is operated for the first time, make sure that the safety valve is working properly. Turn screw (4) of safety valve (1) several rotations to the left until the safety valve releases air. Let the safety valve blow out for only a few seconds. Turn screw (4) to the right until it seats, closing the valve.

The safety valve must never be used for depressurizing the air tank. It could damage the safety valve. The valve is set to the maximum permitted pressure by the manufacturer. Adjustments are not permitted.

Warning! Compressed air can be dangerous. Wear eye protection when blowing air out.

14.3. Replacement of the input filter and prefilter

(Fig.16)

At the lid of the compressors crankcase is an input filter (1) and prefilter (3).

Replacing of the input filter:

- Hand pull the rubber stopper (2).
- Remove used and dirty filter.
- Input new filter and set rubber stopper.

Replacing of the prefilter:

- Hand pull prefilter (3).
- Replace old prefilter with new.


Fig.16

14.4. Replacement of output filter in air dryer

Before beginning, reduce the air pressure in the tank to zero and disconnect the appliance from the mains.

(Fig.17)

In the case of a regular operation of a dryer it is necessary to replace the dryer filter in the upper part or after repair of dysfunction which could pollute them.

- Remove the dryer plug (1) by unscrewing to the left from the dryer head (4).
- Replace filter with the new one (2) and clean the sieve (3).
- It is possible to check or to change the filling
- Refit the dryer plug to the dryer head (4) and tighten it to the right.


Fig.17

14.5. Replacement of filter in condensation and filtration unit

Before beginning, depressurize the air tank to zero and disconnect the appliance from the mains.

(Fig.18)

In the case of a regular operation of a condensation unit it is necessary to replace the filter inside the filter with automatic desludging.

- Release a safety lock (1) on the filter vessel by its pulling downwards, slightly rotate the filter cover (2) to the left and take it out.
- Unscrew the filter holder (3) by its rotation to the left.
- Replace the filter and fix the new one by rotation of the holder to the right back on the filter body.
- Replace the filter cover and secure it by turning to the right until the safety pin locks.

Fig.18

14.6. Replacement of the buoy in the water separator

(Fig.19)

During regular operation of drier is necessary to change the buoy in the water separator.

- Release pressure.
- Dismount separator container.
- Pull out condenser separator
- Release nut of the buoy placed in the bottom part of the container.
- Pull out worn buoy of the separator and replace it for the new one.
- Secure the buoy with the nut in bottom part of the container.
- Place the condenser separator back as shown in the picture.
- Container of the separator insert back and screw in.
- The container is definitely secured in the point of the symbol.


Fig. 19

15. STORAGE

If the compressor will not be used for a prolonged time period, drain any condensate from the air tank. Then turn on the compressor for 10 minutes, keeping the drain valve open (1) (Fig.14, Fig.15). Switch off the compressor by switch (3) at pressure switch (2) (Fig.13), close the drain valve and disconnect the appliance from the mains.

16. DISPOSING OF THE APPLIANCE

- Disconnect the appliance from the mains.
- Release air pressure in the pressure tank by opening the drain valve (1) (Fig.14, Fig.15).
- The components of the product are non-toxic.
- Dispose of the appliance following all environmental regulations.

17. REPAIR SERVICE

Guaranteed and post-guarantee repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

The manufacturer reserves the right to make changes to the appliance without notice. Any changes made will not affect the functional properties of the appliance.

18. SOLVING PROBLEMS

Caution! Before proceeding, depressurize the air tank to zero and disconnect the appliance from the mains.

After the remedy of failure and reverse assembly of drier, it is necessary to perform the regeneration of drier, best at continuous operation of compressor at pressure of 0.6-0.7 MPa for at least 1 hour and then to check air drying.

Troubleshooting can be performed only by qualified personnel.

FAILURE	POSSIBLE CAUSE	REMEDY
Compressor does not start	No voltage in pressure switch	Check voltage in socket Check fuse – replace faulty one Loosen terminal – tighten it Check power cord – replace faulty one Replace motor or re-wind it
	Disconnected winding of motor, damaged thermal protection Faulty capacitor Seizure of piston or another rotary part Pressure switch does not switch on	Replace capacitor Replace damaged parts Check the function of pressure switch
Compressor often switches on	Air leak in pneumatic distribution system Leaking check valve Greater volume of condensed liquid in pressure vessel	Check pneumatic distribution system – seal loose joint Clean valve, replace seals, replace valve Drain condensed liquid
Prolonged running of compressor	Air leak in pneumatic distribution system Worn piston ring Contaminated input filter and prefilter Dirty filter in the dryer Defective solenoid valve	Check pneumatic distribution system – seal loose joint Replace worn piston ring Replace contaminated filters with the new ones Change the outlet filter and inspect dessicant Repair or change the valve
Compressor is noisy (knocking, metal noises)	Damaged bearing of piston, piston rod, motor bearing Loose or cracked spring	Replace damaged bearing Replace damaged spring
Dryer doesn't dry (condensed water in the tank)	Devaluated or very polluted dessicant	Change the dessicant and filters
	low operating pressure	reduce the air offtake, check the source efficiency, check possible leaks in the distribution
	inoperative solenoid valve	repair or change the valve
	clogged nozzle of the regenerating air	clean or replace the nozzle, use correct size of the nozzle (see Product maintenance)
	inoperative cooler ventilator	replace ventilator check supply of electric energy
Drying unit is noisy or sounds annoying	leak of white liquid through the solenoid valve	dismantle the chamber, replace drying substance, bottom filter and sealing and inspect tightness, apply soapy water onto O-rings of the nuts
	defective solenoid valve	replace the valve
	damaged damping substance in the condensate vessel	replace damping substance or the vessel
	damaged pressure hose	replace pressure hose

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