**DK50 2V/50**  
**DK50 4VR/50**

INSTALLATION, OPERATION AND MAINTENANCE MANUAL  
DIE ANLEITUNG ZUR INSTALLATION, BEDIENUNG UND WARTUNG  
ИНСТРУКЦИЯ ПО УСТАНОВКЕ, ОБСЛУЖИВАНИЮ И УХОДУ  
NÁVOD NA INŠTALÁCIU, OBSLUHU A ÚDRŽBU

Ⓔ

Ⓓ

Ⓡ

Ⓚ

**ekom**<sup>®</sup>

EN	<a href="#"><u>CONTENTS..... 3</u></a>
DE	<a href="#"><u>INHALT..... 20</u></a>
RU	<a href="#"><u>СОДЕРЖАНИЕ..... 38</u></a>
SK	<a href="#"><u>OBSAH..... 57</u></a>

**CONTENTS**

**IMPORTANT INFORMATION..... 4**

- 1. CE MARKING..... 4
- 2. WARNINGS..... 4
- 3. ALERT NOTICES AND SYMBOLS..... 5
- 4. STORAGE AND TRANSPORT ..... 5
- 5. TECHNICAL DATA..... 6
- 6. PRODUCT DESCRIPTION ..... 7
- 7. FUNCTION ..... 7

**INSTALLATION ..... 10**

- 8. USE ..... 10
- 9. INSTALLATION ..... 10
- 10. WIRING DIAGRAMS ..... 12
- 11. PNEUMATIC DIAGRAM..... 15
- 12. FIRST OPERATION ..... 16

**OPERATION ..... 16**

- 13. SWITCHING THE COMPRESSOR ON ..... 16

**MAINTENANCE ..... 17**

- 14. MAINTENANCE SCHEDULE..... 17
- 15. MAINTENANCE ..... 17
- 16. STORAGE ..... 19
- 17. DISPOSING OF THE APPLIANCE ..... 19
- 18. REPAIR SERVICE ..... 20
- 19. SOLVING PROBLEMS..... 20

**PARTS LIST..... 79**

## **IMPORTANT INFORMATION**

### **1. CE MARKING**

Products labeled with the CE mark of compliance meet the safety guidelines (93/42/EEC) of the European Union.

### **2. WARNINGS**

#### **2.1. General warnings**

- This Installation, Operation and Maintenance Manual is a part of the appliance and must be kept with the compressor. Careful review of this manual will provide the information necessary for correct operation of the appliance.
- The safety of operating personnel and trouble-free operation of the appliance are guaranteed only if original parts are used. Only accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer can be used.
- If any other accessories or consumable materials are used, the manufacturer cannot be held responsible for the safe operation of the appliance. This guarantee does not cover damages originating from the use of accessories or consumable material other than those specified or suggested by the manufacturer.
- The manufacturer guarantees the safety, reliability and function of the appliance only if:
  - Installation, new settings, amendments, extensions and repairs are performed by the manufacturer or its representative, or a service provider authorized by the manufacturer
  - The appliance is used in accordance with this Installation, Operation and Maintenance Manual
- The manufacturer reserves all rights for the protection of its wiring diagrams, methods and names.
- Translation of Manual for Installation, Operation and Maintenance is carried out in accordance with the best knowledge. In the case of ambiguities, the Slovak version of the text prevails.

#### **2.2. General safety warnings**

The manufacturer developed and designed the equipment in such a way so that any risks were excluded if it is used according to intention. The manufacturer considers it to be its obligation to describe the following safety measures in order to exclude residual damages.


- Operation of the appliance must be in compliance with all local codes and regulations.
- Original packaging should be kept for the return of the appliance. Only the original packaging ensures protection of the appliance during transport. If it is necessary to return the appliance during the guarantee period, the manufacturer is not liable for damages caused by improper packaging.
- Each time the appliance is used, the operator must make sure that it is functioning correctly and safely.
- The user must fully understand the operation of the appliance.
- The product is not intended for operation in areas with a risk of explosion.
- If any problem occurs during use of the appliance, the user must inform his supplier immediately.

#### **2.3. Electrical system safety warnings**

- The appliance must be connected to earth (grounded).
- Before the appliance is plugged in, make sure that the mains voltage and mains frequency stated on the appliance are the same as the power mains.
- Prior to putting into operation it is necessary to check for possible damage of the equipment and connected air and electric distributions. Damaged pneumatic and electric lines must be immediately replaced.
- Immediately disconnect the appliance from the mains (pull out mains plug) if a technical failure occurs.
- During repairs and maintenance, ensure that:
  - The mains plug is pulled out from the socket
  - Pressure pipes are vented and pressure is released from the air tank.
- The appliance must be installed by an approved, qualified technician.

### 3. ALERT NOTICES AND SYMBOLS

In the Installation, Operation and Maintenance Manual and on packaging and product, the following labels or symbols are used for important information:

	Attention, see instructions for use
	Caution, risk of electric shock
	Consult instructions for use
	CE mark of compliance
	Compressor is remote-controlled and may start without warning
	Caution! Hot surface
	Earth (ground) connection
	Terminal for ground connection
	Fuse
	Alternating current
	Handling mark on package – Fragile, handle with care
	Handling mark on package – This way up (vertical position of cargo)
	Handling mark on package – Protect against moisture
	Handling mark on package – Temperature during storage and transport
	Handling mark on package – Limited stacking
	Mark on package – Recyclable material

### 4. STORAGE AND TRANSPORT

The compressor is shipped in cardboard that protects the appliance from damage during transport.

**Caution! For transport, always use the original packaging and secure the compressor in the upright position.**

Protect the compressor from humidity and extreme temperatures during transport and storage. A compressor in its original packaging can be stored in a warm, dry and dust-free area. Do not store near any chemical substances.

Keep packaging material if possible. If not, please dispose of the packaging material in an environmentally friendly way and recycle if possible.

**Caution! Before moving or transporting the compressor, release all the air pressure from the tank and hoses and drain the condensed water.**

**5. TECHNICAL DATA**

	DK50 2V/50	DK50 4VR/50	DK50 2V/50S	DK50 4VR/50S
Nominal voltage / frequency V / Hz	230 / 50 3x400/50	3x400/50	230 / 50 3x400/50	3x400/50
Efficiency of compressor at over-pressure 6 bar Lit.min <sup>-1</sup>	134	270	134	270
Efficiency of compressor with dryer at over-pressure 6 bar Lit.min <sup>-1</sup>	95(**)	200(**)	95(**)	200(**)
Maximal current A	7.4 4.5	4.9	7.6 4.6	5.3
Maximal current of compressor with dryer A	7.6 4.7	5.2	7.8 4.7	5.6
Motor performance kW	1.1 1.5	2.2	1.1 1.5	2.2
Air tank capacity Lit.	50	50	50	50
Pressure range bar	6,0 – 8,0	6,0 – 8,0	6,0 – 8,0	6,0 – 8,0
Maximum operating pressure of safety valve bar	12,0	12,0	12,0	12,0
Sound level L <sub>pfA</sub> [dB]	71	75	56	59
Mode of operation of compressor	CONTINUAL S 1	CONTINUAL S 1	CONTINUAL S 1	CONTINUAL S 1
Mode of operation of compressor with dryer	CONTINUAL S 1	CONTINUAL S 1	CONTINUAL S 1	CONTINUAL S 1
Dimensions of compressor / of compressor with dryer w x l x h mm	595x475x770 / 595x585x770	595x475x780 / 595x585x780	750x715x1015	750x715x1015
Dimensions of compressor packaged in cardboard w x l x h mm			945x890x1205	945x890x1205
Weight of compressor / of compressor with dryer kg	56 / 70	70 / 90	108 / 123	120 / 140
Weight of compressor / of compressor with dryer Packaget in cardboard kg				
Drying point of compressor Atmospheric condensation point	-20°C			
Version EN 60 601-1	Appliance of type B, class I			





(\*\*) Part of the efficiency of compressor (cca-20%) is used for regeneration of the dryer  
 - Air outgoing from dryer M3, M2 is filtered using 5µm filter

**Climatic conditions during storage and transport**

Temperature : -25°C to +55°C, 24 h to +70°C  
 Relative air humidity : 10% to 90 % (no condensation)

**Climatic operation conditions**

Temperature : +5°C to +40°C  
 Relative air humidity : 70%

## 6. PRODUCT DESCRIPTION

### 6.1. Model variations and their uses

Compressors are the source of clean, oil-free compressed air used to drive dental appliances and equipment.

Compressors models are designed for the following uses:

**Dental compressors DK50 2V/50** - are designed for independent placement of the compressor in any area.

**Dental compressors DK50 2V/50/M** - are designed for independent placement of the compressor in any area and feature an adsorption air dryer M3.

**Dental compressors DK50 2V/50S** - feature soundproof boxes suitable for placing in the dentist's surgery.

**Dental compressors DK50 2V/50S/M** - feature soundproof boxes and an adsorption air dryer M3.

**Dental compressors DK50 4VR/50** - These compressors can be situated in areas where their operation and noise will not disturb their immediate surroundings. They are capable of providing compressed air for several dentists' surgeries.

**Dental compressors DK50 4VR/50/M** - are equipped with an adsorption air dryer M2.

**Dental compressors DK50 4VR/50S** - are placed in compact soundproof boxes and are capable of providing compressed air for several dental units.

**Dental compressors DK50 4VR/50S/M** - are placed in compact soundproof boxes and are equipped with an adsorption air dryer M2.

**Without additional filtration equipment, the compressed air from a compressor is not suitable for the operation of breathing appliances or similar equipment.**

## 7. FUNCTION

### Compressor (Fig.1)

The compressor (1) draws in air through a filter (8) and compresses it through a check valve (3) into an air tank (2). The connected apparatus draws the compressed air from the air tank until the pressure drops to a default preset level on the air-pressure switch (4) switching the compressor on. The compressor again compresses air into the nozzle until the maximum pressure is reached and the compressor switches off. After compressor aggregate is switched off, pressure hose shall be pressure-release solenoid valve (21). Safety valve (5) prevents the pressure in air chamber from rising above the maximal allowed value. The drain valve (7) releases the condensate from the air nozzle. Compressed, clean air free from oil traces is stored in the air tank ready for use.

### Compressor with dryer M2, M3 (Fig.2, Fig.3)

Compressor aggregate (1) sucks the atmospheric air through input filter (8) and it supplies it compressed to the cooler. Air further proceeds via condensate trap through control valve (12), through drier chamber (9) with adsorber (18), where humidity is captured, through output filter (16) and a check valve (3) to air chamber (2). After each switching off of the compressor using the pressure switch the adsorber is regenerated while the drier chamber is depressurized. Air escapes from adsorption chamber through the open solenoid valve (10), while it is concurrently blown with dried air. Air drying takes place in one chamber and regeneration takes place in the other chamber. The mode of the chambers is changed in regular cycles

and the drying and regeneration process takes place in the chambers in a reversal order. The compressed, dry and clean air free from oil traces and is ready in the air chamber for further use.

**Compressor box (Fig.4)**

The soundproof box is compact yet allows sufficient exchange of cooling air. It can be placed in a dentist's office. The ventilator (13) under the aggregate of a compressor provides cooling of compressor and it is in operation at the same time with an engine of the compressor. After a prolonged operation of the compressor, when temperature in a casing rises above 40°C, the cooling ventilator of the casing (22) shall be switched on automatically. After the area in housing is cooled down under ca 32°C, the ventilators shall be automatically switched off.


**Make sure that nothing impedes the free flow of air under and around the compressor. Never cover the hot air outlet on the top back side of the case.**

**If placing the compressor on a soft floor such as carpet, create space for ventilation between the base and floor or the box and floor, e.g. underpin the footings with hard pads.**

**Fig.1 - Compressor**

1. Compressor motor
2. Air tank
3. Check valve
4. Pressure switch
5. Safety valve
6. Manometer
7. Condenser outlet
8. Input filter
9. Dryer chamber
10. Output solenoid valve of dryer
11. Cooler
12. Control valve
13. Compressor fan
14. Stopper
15. O-ring
16. Output filter
17. Sieve
18. Dessicant
19. Water separator
20. Bottle
21. Solenoid valve
22. Box fan (DK50 4VR - 2x)
23. Thermo switch
24. Switch
25. Manometer
26. Konektor skrinky
27. Hose of manometer
28. Magnetic bottle holder
29. Connecting reinforcement
30. Pulling system
31. Wheels
32. Wall stopper
33. Check valve

**Fig.2 – Compressor with dryer MONZUN – M3**


Fig.3 – Compressor with dryer MONZUN – M2



Fig.4 - Box

## INSTALLATION

### 8. USE

- The appliance must be installed and operated in a dry, well ventilated and dust-free area where ambient temperature is within the range of +5°C to +40°C and relative air humidity does not exceed 70%. Otherwise, failure-free operation of the compressor cannot be guaranteed. The compressor must be installed so that it is accessible at all times for operation and maintenance. Please ensure that the appliance label is accessible.
- The appliance must stand on a flat, sufficiently stable base. See paragraph 5 (Technical data) when positioning or lifting the compressor.
- Compressors cannot be exposed to outdoor environments. The appliance cannot be used in moist or wet environments. Do not use the compressor in the presence of explosive gases, dust or combustible liquids.
- Before connecting the compressor to medical equipment, the supplier must confirm that it meets all requirements for its use. Refer to the technical data of the product for this purpose. When a unit is to be built-in, classification and evaluation of compatibility must be done by the manufacturer or supplier of the product to be used.
- Any use other than that described in this manual is not covered by the guarantee, and the manufacturer is not liable for any damages that may result. The operator/user assumes all risk.

### 9. INSTALLATION

Only qualified personnel can install and start up the appliance and train operating personnel in its correct use and maintenance. Installation and training of all operators shall be confirmed by the installer's signature on the certificate of installation.

Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.


Caution! When in operation, the compressor is hot. Burns or fire may result if contact is made by the operator or any flammable material.

Electric cord for connection to electric mains and air hoses may not be broken. The power cord may not be exposed to pulling, pressure and excessive heat.

#### 9.1. Placement of the compressor

Handling


Fig.5 Unpacking

**Dental compressor DK50 2V/50, DK50 2V/50/M, DK50 4VR/50, DK50 4VR/50/M (Fig.5)**

After removing all packaging material, place the product on the floor and remove stabilization parts Y (Detail A). Connect output hose with end-piece to the appliance. Plug the mains cord plug into a socket.

**Dental compressor in box DK50 2V/50S, DK50 4VR/50S (Fig.4, Fig.5)**

After unpacking, place the product onto the floor in a room, release it from packaging materials and remove fixation parts (Y) - detail A. Place a wall end stop (32) 2 pcs on compressor casing in the rear, upper part of the casing and fit in the casing on the required place. The end-stops shall ensure the sufficient distance of the casing from the wall for thorough ventilation. For fitting the compressor in the place to the casing it is necessary to disassemble the casing door and take off the connecting reinforcement (29) in the front bottom part of the casing. Protrude the pressure hose under the casing and fix it to the appliance in a suitable manner. Grasp the compressor at its handle and using the transport mechanism (30) and the built-in castors (31) place it to the casing. Embed the hose (27) of manometer (25) of the casing to fast-on coupling on the compressor, place back the connecting reinforcement (29) and connect output pressure hose to the compressor.

Connect the connector (26) of casing to the compressor and connect the electric mains plug to mains socket.

**In case of disassembling the compressor it is necessary to disconnect the connector of soundproofing box by using the screwdriver.**

(Fig.6)

Fig.6

**Dental compressor in DK50 2V/50S/M, DK50 4VR/50S/M (Fig.4, Fig.5)**

After removing all packaging material, place the product on the floor and remove stabilization parts Y (Detail A). Place the compressor into housing similarly as in the previous paragraph. Prior to placing the compressor into housing, protrude house for condensate drain via hole in housing and connect it to a bottle (20). Magnetic holder (28) with a vessel (20), for entrapping the condensate from the dryer may be fitted onto any vertical part of casing, or from front on its door. When fixing the holder with a vessel at the housing side it is necessary to consider a space of at least 11 cm between the housing and furniture. Distance smaller than the specified one may cause problem with handling of the vessel.

**9.2. Compressed air outlet**

(Fig.7)

Lead the pressure hose from the output of compressed air (1) to the appliance – dental set.


Fig.7

**9.3. Electrical connection**

**Plug the electrical cord into the mains.**

**The appliance is equipped with a grounded plug. Make sure this connection complies with local electrical codes. The mains voltage and frequency must comply with the data stated on the appliance label.**

(Fig.8)

- Keep the socket easily accessible to ensure that in an emergency the appliance can be safely disconnected from the mains.
- Connection to the power distribution box must be max.16 A.
- The connection of the earth ground pin  $\varnothing$  6mm (1) with other appliances must be completed in accordance with local electrical codes. The female socket (2), which is not included in the standard set, is an optional accessory.

Fig.8


**Electrical cable may not contact the hot parts of a compressor. Insulation could be damaged!**

**If any electrical cord or air hose is damaged it must be replaced immediately.**

10. WIRING DIAGRAMS


3N/PE - 400/230 V 50 Hz  
MAINS TN-S (TN-C-S)  
ELECTRIC OBJECT 1ST CAT.  
TYPE B

- SP Pressure switch
- EV1, EV2 Fan of compressor
- ST1 Thermo switch
- EV3 Fan of dryer
- M1 Motor of compressor
- FA1 Breaker
- YV1 Solenoid valve
- YV2, YV4 Solenoid valve of dryer - OUT
- KM1 Contactor
- YV3 Control valve
- X1 Terminal
- KT1 Time relay
- PH1\* Hour counter
- XC1 Connector
- Q Switch
- HL1 Glowlamp
- ST2 Thermoswitch of box
- EV4, EV5 Fan of box


DK50 4VRS / 50

3N/PE - 400/230 V 50 Hz  
MAINS TN-S (TN-C-S)  
ELECTRIC OBJECT 1ST CAT.  
TYPE B


DK50 4VRS / 50 / M2

3N/PE - 400/230 V 50 Hz  
MAINS TN-S (TN-C-S)  
ELECTRIC OBJECT 1ST CAT.  
TYPE B


DK50 4VR / 50

3N/PE - 400/230 V 50 Hz  
MAINS TN-S (TN-C-S)  
ELECTRIC OBJECT 1ST CAT.  
TYPE B


DK50 4VR / 50 / M2


30. 11. 2012

DK50 2V / 50


DK50 2VS / 50

DK50 2V / 50 / M3


- SP Pressure switch
- EV1 Fan of compressor
- ST1 Thermo switch
- EV3 Fan of dryer
- M1 Motor of compressor
- FA1 Breaker
- YV1, YV2, YV4 Solenoid valve
- YV3 Solenoid valve of dryer - OUT
- KM1 Contactor
- YV3 Control valve
- X1 Terminal
- KT1 Time relay
- PH1\* Hour counter
- XC1 Connector
- Q Switch
- HL1 Glowlamp
- ST2 Thermo switch of box
- EV4 Fan of box
- CB1 Capacitor

DK50 2VS / 50 / M3


3/N/PE ~ 400/230 V 50 Hz  
MAINS TN-S [TN-C-S]  
ELECTRIC OBJECT 1ST CAT.  
TYPE B

\* - only for 10 bar


3/N/PE ~ 400/230 V 50 Hz  
MAINS TN-S [TN-C-S]  
ELECTRIC OBJECT 1ST CAT.  
TYPE B

\* - only for 10 bar

3/N/PE ~ 400/230 V 50 Hz  
MAINS TN-S [TN-C-S]  
ELECTRIC OBJECT 1ST CAT.  
TYPE B


TIME RELAY  
KT1 - 90 sek.

\* - only for 10 bar

- SP Pressure switch
- EV1 Fan of compressor
- ST1 Thermo switch
- EV3 Fan of dryer
- M1 Motor of compressor
- FA1 Breaker
- YV1 Solenoid valve
- YV2 , YV4 Solenoid valve of dryer - OUT
- KM1 Contactor
- YV3 Control valve
- X1 Terminal
- KT1 Time relay
- PH1\* Hour counter
- XC1 Connector
- Q Switch
- HL1 Glowlamp
- ST2 Thermo switch of box
- EV4 Fan of box

3/N/PE ~ 400/230 V 50 Hz  
MAINS TN-S [TN-C-S]  
ELECTRIC OBJECT 1ST CAT.  
TYPE B


TIME RELAY  
KT1 - 90 sek.


\* - only for 10 bar

DK50 2VS / 50 / M3

**11. PNEUMATIC DIAGRAM**

DK50 2V/50/M

DK50 2V/50

- 1 Input filter (for 4V/50- 2x)
- 2 Compressor motor
- 3 Compressor Fan
- 4 Solenoid valve
- 5 Safety valve
- 6 Check valve
- 7 Pressure switch
- 8 Air tank
- 9 Manometer
- 10 Output solenoid valve of dryer
- 11 Cooler
- 12 Cooler Fan
- 13 water separator
- 14 Relief valve
- 15 control valve
- 16 Dryer
- 17 Jet

## 12. FIRST OPERATION

(Fig.98)

- Make sure that all stabilizers used during transport were removed.
- Check that all pressurized air line connections are secure.
- Connect to the mains.
- Start compressor at pressure switch (2) by turning switch (3) to position "I."
- For kompressor in the box turn the switch (24) (Fig.4) at the front part of the soundproof box to the position "I" – green light indicates that the appliance is on.

**Compressor** - At first operation the air tank is pressurized until it reaches a preset level when the compressor automatically switches off. As the air is used, the compressor works in automatic mode, switched on or off by the pressure switch.

**Compressor with dryer** - In addition, adsorption dryer takes out humidity from the previous compressed air during operation in the device and it blows out the entrapped condensate via condensate outlet on the dryer, which is audible as a short hiss when compressor stops, or during the operation while switching the dryer chambers.

**The compressor is not equipped with an emergency power supply.**

### OPERATION

**In case of emergency, disconnect the compressor from the mains (pull out the mains plug).**

**The compressor has hot surfaces.  
Burns or fire may result if contact is made.**

**During prolonged operation of the compressor, the temperature in the box may increase to over 40°C. At this point the cooling fan automatically switches on. After cooling the space to under 32°C, the ventilator switches off.**

**Automatic start: when pressure in the tank drops to the pressure switch's lower limit level, the compressor automatically switches on. The compressor automatically switches off after reaching the pressure switch's upper limit level.**


### **Compressor with dryer**

- It is forbidden to alter the working pressures of pressure switch set by manufacturer. The operation of the compressor at working pressure lower than the switching pressure demonstrates the overload of the compressor (high air consumption) by the appliance, leakages in pneumatic distributions, failure of aggregate or drier.
- Prior connecting drier to air chamber, that was used with compressor without drier, it is necessary to clean interior surface of air chamber and perfectly remove condensed liquid. Then interconnect electric part of drier with compressor according to wiring diagram in accord with valid regional regulations.

## 13. SWITCHING THE COMPRESSOR ON

(Fig.9)

Switch on the compressor at the pressure switch (2) by turning the knob (3) to position "I.", (for compressor in the box switch (24) Fig.4 on the front part of the compressor box), The compressor sends pressurized air to the air tank. As the compressed air is used, the pressure in the air nozzle drops to a preset level, the compressor switches on and the air nozzle fills with compressed air. After reaching the cutoff pressure the compressor turns off automatically and the cycle is repeated. Check the value of switching-on and switching-off pressure on pressure gauge. The values may be within a tolerance of  $\pm 10\%$ . Air pressure in air chamber must not exceed maximal permitted operation pressure. Fig.9

**Never tamper with the pressure switch (2). Adjustments are not allowed. The pressure switch (2) has been set by the manufacturer and further setting of switching on and off pressure may be carried out only by a qualified expert trained by the manufacturer.**

Size of regeneration nozzle for dryer M2, M3

Compressor	Jet size	Closing pressure - cut-off pressure	Amount of dried air	Amount of regenerated air
DK50 2V/50/M3	Ø 0,6 mm	6 - 8 bar	140 l/min.	30 l/min./ 6bar
DK50 4VR/50/M2	Ø 0,85 mm	5 - 8 bar	280 l/min.	50 l/min./ 5bar

**MAINTENANCE**

**14. MAINTENANCE SCHEDULE**

**Notice!**

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

Maintenance that must be performed	Chapter	Time interval	Performed by
<ul style="list-style-type: none"> <li>• Release condensate Compressor without air drier At high air humidity Compressors with air drier</li> </ul>	15.1	1 x week 1 x day 1 x week, check function	operating staff operating staff operating staff
<ul style="list-style-type: none"> <li>• Check safety valve</li> </ul>	15.2	1 x year	qualified technician
<ul style="list-style-type: none"> <li>• Replacement of the input filter and prefilter</li> </ul>	15.3	1 x 2 years or after 5000 hours	qualified technician
<ul style="list-style-type: none"> <li>• Replacement of filter in dryer</li> </ul>	15.4	1 x year	qualified technician
<ul style="list-style-type: none"> <li>• Replacement of the buoy in the water separator</li> </ul>	15.5	1 x year	qualified technician
<ul style="list-style-type: none"> <li>• Check tightness of joints Overall examination of device</li> </ul>	Service documentation	1 x year	qualified technician
<ul style="list-style-type: none"> <li>• Perform "Repeated Test" according to EN 62353</li> </ul>	14	1 x 2 years	qualified technician

**15. MAINTENANCE**

Repair work beyond normal maintenance can be performed only by qualified personnel or the manufacturer's representative.  
Use only spare parts and accessories approved by the manufacturer.

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

For permanently high efficiency of drying, it is necessary to maintain the whole appliance, and mainly ventilator clean – regularly clean the surface of ventilator and cooling fins of cooler.

TO ENSURE THAT THE COMPRESSOR WORKS CORRECTLY, PERFORM THE FOLLOWING MAINTENANCE TASKS AT REGULAR INTERVALS (CHAPTER 14):.

### 15.1. Condensation drain valve

#### Compressors (Fig.10)

During regular use, release condensation from the pressure tank. Switch off the compressor at the mains. Reduce air pressure in the appliance to max. 1 bar by releasing air via a connected device. Place the vessel under release valve and open the drain valve (1). Wait until condensation is fully drained from the pressure tank. Close drain valve (1).


Fig.10

#### Compressors with air dryer

In the case of a regular operation condensate is automatically excreted via air dryer and it is entrapped in a bottle. Take out the bottle from a holder, release cap and pour out the condensate.

### 15.2. Safety valve check

(Fig.8)

When the compressor is operated for the first time, make sure that the safety valve is working properly. Turn screw (4) of safety valve (1) several rotations to the left until the safety valve releases air. Let the safety valve blow out for only a few seconds. Turn screw (4) to the right until it seats, closing the valve.

**The safety valve must never be used for depressurizing the air tank. It could damage the safety valve. The valve is set to the maximum permitted pressure by the manufacturer. Adjustments are not permitted.**

**Warning! Compressed air can be dangerous. Wear eye protection when blowing air out.**

### 15.3. Replacement of the input filter and prefilter

(Fig.11)

At the lid of the compressors crankcase is an input filter (1) and prefilter (3).

Replacing of the input filter:

- Hand pull the rubber stopper (2).
- Remove used and dirty filter.
- Input new filter and set rubber stopper.

Replacing of the prefilter:

- Hand pull prefilter (3).
- Replace old prefilter with new.


Fig.11

### 15.4. Replacement of output filter in air dryer

**Before beginning, reduce the air pressure in the tank to zero and disconnect the appliance from the mains.**

(Fig.12)

In the case of a regular operation of a dryer it is necessary to replace the dryer filter in the upper part.

- Remove the dryer plug (1) by unscrewing to the left from the dryer head (4).
- Put in a new filter (2).
- Refit the dryer plug to the dryer head (4) and tighten it to the right.


Fig.12

**15.5. Replacement of the buoy in the water separator**

(Fig.13)

During regular operation of drier is necessary to change the buoy in the water separator.

- A) Release pressure.
- B) Dismount separator container.
- C) Pull out condenser separator
- D) Release nut of the buoy placed in the bottom part of the container.
- E) Pull out worn buoy of the separator and replace it for the new one.
- F) Secure the buoy with the nut in bottom part of the container.
- G) Place the condenser separator back as shown in the picture.
- H) Container of the separator insert back and screw in.
- I) The container is definitely secured in the point of the symbol.


Fig.13

Water separator	Order number	Buoy	Order number
WS010BBFX	025200119	EF1	025200146

**16. STORAGE**

If the compressor will not be used for a prolonged time period, drain any condensate from the air tank. Then turn on the compressor for 10 minutes, keeping the drain valve open (1) (Fig.9). Switch off the compressor by switch (3) at pressure switch (2) (Fig.8), close the drain valve and disconnect the appliance from the mains.

**17. DISPOSING OF THE APPLIANCE**

- Disconnect the appliance from the mains.
- Release air pressure in the pressure tank by opening the drain valve (1) (Fig.9).
- The components of the product are non-toxic.
- Dispose of the appliance following all environmental regulations.

## 18. REPAIR SERVICE

Guaranteed and post-guarantee repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

The manufacturer reserves the right to make changes to the appliance without notice. Any changes made will not affect the functional properties of the appliance.

## 19. SOLVING PROBLEMS

**Caution! Before proceeding, depressurize the air tank to zero and disconnect the appliance from the mains.**

**After the remedy of failure and reverse assembly of drier, it is necessary to perform the regeneration of drier, best at continuous operation of compressor at pressure of 0.6-0.7 MPa for at least 1 hour and then to check air drying.**

Troubleshooting can be performed only by qualified personnel.

FAILURE	POSSIBLE CAUSE	REMEDY
Compressor does not start	No voltage in pressure switch  Disconnected winding of motor, damaged thermal protection Faulty capacitor Seizure of piston or another rotary part Pressure switch does not switch on	Check voltage in socket Check fuse – replace faulty one Loosen terminal – tighten it Check power cord – replace faulty one Replace motor or re-wind it  Replace capacitor Replace damaged parts Check the function of pressure switch
Compressor often switches on	Air leak in pneumatic distribution system Leaking check valve Greater volume of condensed liquid in pressure vessel	Check pneumatic distribution system – seal loose joint Clean valve, replace seals, replace valve Drain condensed liquid
Prolonged running of compressor	Air leak in pneumatic distribution system Worn piston ring Contaminated input filter and prefilter  Dirty filter in the dryer  Defective solenoid valve	Check pneumatic distribution system – seal loose joint Replace worn piston ring Replace contaminated filters with the new ones Change the outlet filter and inspect dessicant Repair or change the valve
Compressor is noisy (knocking, metal noises)	Damaged bearing of piston, piston rod, motor bearing Loose or cracked spring	Replace damaged bearing  Replace damaged spring
Dryer doesn't dry (condensed water in the tank)	Devaluated or very polluted dessicant	Change the dessicant and filters
	low operating pressure	reduce the air offtake, check the source efficiency, check possible leaks in the distribution
	inoperative solenoid valve	repair or change the valve
	clogged nozzle of the regenerating air	clean or replace the nozzle, use correct size of the nozzle (see Product maintenance)
	inoperative cooler ventilator	replace ventilator check supply of electric energy
Drying unit is noisy or sounds annoying	leak of white liquid through the solenoid valve	dismantle the chamber, replace drying substance, bottom filter and sealing and inspect tightness, apply soapy water onto O-rings of the nuts
	defective solenoid valve	replace the valve
	damaged damping substance in the condensate vessel	replace damping substance or the vessel
	damaged pressure hose	replace pressure hose

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