

CHEESE UNI - CHIRANA's news designed for unique solutions



New design intended for work with lying patient - in position "12". CHEESE UNI can be equipped from a minimal version suitable for orthodontic rooms up to standard version of treatment unit. CHIRANA CHEESE UNI is designed for left and right handed dentists with rotary dentist and assistant table to right or to left depending on dentist's need.















Dentist table

5 places: syringe/ turbine/ micromotor/ scaler / big suction/ small suction



Control panel

- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control (by foot control or manual from keyboard)
- Disinfection programm activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control



... with assistant arm

you can place big and small suction hoses on the assistan table and save 2 places on the dentist table

... Or with tray table

1x silicone or 2x silicone tray table behind the chair





Optional suction systems



Suction systems:

- DURR MSBV valve
- Ready to connect to central wet suction
- Ejector suction

Hygiene:

- Clean water bottle system
- Disinfection system
- Clean water bottle + disinfection system
- DCI bottles
- Easy to clean Durr filter
- Easy to clean a suction filter
- Easy to clean an oil-catcher













FARO ALYA LED

FARO MAIA LED

adjustable light intensity 3,000 - 35,000 lux

LED - A adjustable light intensity 6,000 - 38,000 lux

FARO EDI

adjustable light intensity 8,000 - 25,000 lux



adjustable light intensity 3,000 - 50,000 lux









Multifunctional foot control FCR1-7- WL Multifunctional foot control FCR1-7

Multifunctional foot pedal

Foot control with spray

Foot switch (only istruments activation)



